



- ① Ergonomic Ear Hook
- ② Microphone
- ③ Volume Control
- ④ Battery compartment (Type 13)
- ⑤ Acoustic Filter

DESCRIPTION

Behind-The-Ear air conduction hearing aid (BTE), Battery 13, programmable through fitting software, WDRC strategy.

INTENDED USE

The Microson m2 BTE is indicated to compensate from moderate to severe hearing loss (transmissive, mixed or sensorineural). It is not intended for children. See Fitting Range⁽¹⁾


Features

- ✓ Matrix 131/58@ 2cc
- ✓ 100% Digital Technology
- ✓ Digitally Programmable
- ✓ 2 Programs
- ✓ 6 WDRC Independent Channels **New!**
- ✓ Autophone
- ✓ Telecoil
- ✓ Digital Volume Control
- ✓ 6 MPO Control Channels
- ✓ 12 Frequency Bands
- ✓ Automatic Feedback Canceller
- ✓ Nanoproof Coating⁽²⁾
- ✓ Automatic Noise Reduction up to -12dB
- ✓ Power on Delay
- ✓ Low Battery Indicator
- ✓ Battery 13 Type - PR48 (EN 60086)
- ✓ Small size
- ✓ Suitable for mobile phones⁽³⁾

Requirements


- 89600, Fitting Software Microson CODA e-STUDIO 6 (6.5.0 or higher)
- 53781, 4 Pin Hi-Pro Cable Right
- 53832, 4 Pin Hi-Pro Cable Left
- 79345, Battery 13 Flex Adapter
- 66183, NOAHLINK^A Programming Interface (Kernel v. 1.55.03)
- 73194, HI-PRO^B USB Hearing Instrument Programmer (Firmware 3.00 or higher)
- 88616, HI-PRO^{B2} Hearing Instrument Programmer (Firmware 4.00 or higher)

Weight & Dimensions



Weight* excluding battery: 2,31 gr.
Weight* including battery: 3,09 gr
Dimensions in millimeters (mm)
*Including Damped Ear Hook

(1) Fitting Range



⁽²⁾Nanometric Coating protection against intrusion of particles
⁽³⁾In compliance with IEC 60118-13

DS-0040-002-EN
Rev.C 2017/01/27



m2 BTE
1 of 5

^(A)NOAH & NOAHLINK are licensed products and registered trademarks of HIMSA A/S in Denmark.
^(B)HI-PRO is a registered trademark of GN Otometrics A/S in Denmark.

	Acoustic Data	IEC 60118-7:2005	IEC 60118-0:1993/ A1:1994
OUTPUT	OSPL ¹ 90 Peak (dB _{SPL})	131	136
	OSPL90 Peak Frequency (Hz)	1400	1500
	HFA ² -OSPL90 / RTF ³ -OSPL90 (dB _{SPL})	123	134
GAIN	HFA-FOG ⁴ (dB)	49	57
	RTF-FOG (dB)	50	59
	FOG (dB)	58	61
	FOG Frequency (Hz)	1400	1500
	RTG ⁵ (dB)	45	57
NOISE	Equivalent Input Noise (dB _{SPL})	21	19
AGC ⁶	Attack Time (ms)	1	1
	Release Time (ms)	8	6
TELECOIL	HFA-SPLI ⁷ / RTF-SPLI (dB _{SPL})	105	118
	HFA-FOG-MASL ⁸ / RTF-FOG-MASL @ 1 mA/m (dB _{SPL})	79	88
	500 Hz @ 56 mA/m (THD %)	2.1	4.4
	800 Hz @ 56 mA/m (THD %)	5.9	N/A
	1600 Hz @ 56 mA/m (THD %)	0.8	3.1
DISTORTION	500 Hz @ 70 dB _{SPL} (THD %)	1.5	3.6
	800 Hz @ 70 dB _{SPL} (THD %)	0.9	1.9
	1600 Hz @ 65 / 70 dB _{SPL} (THD %)	0.2	0.4
CONSUMPTION	Current Drain (mA)	0.90	0.88
FREQUENCY LIMITS	f ₁ (Hz)	<100	200*
	f ₂ (Hz)	6500	6700*
		IEC 60318-5:2006	IEC 60318-4:2010

¹OSPL: Output Sound Pressure Level

²HFA: High Frequency Average

³RTF: Reference Test Frequency (1600 Hz)

⁴FOG: Full On Gain

⁵RTG: Reference Test Gain

⁶AGC: Automatic Gain Control

⁷SPLI: Sound Pressure Level Inductive

⁸MASL: Magneto Acoustical Sensitivity Level

*According to DIN 45605 standard


Power Source: 1.3V Battery Simulator

Measurements with an undamped ear hook


Configurable parameters

	Parameter	Minimum	Maximum	Step	Default Parameters	Units
WDRC CHANNELS 1 to 6	Low Level Gain	0	20	2	20	dB
	High Level Gain	0	26	2	26	dB
	Compression TK	45	69	2	55	dB _{SPL}
	Compression Ratio	1:1	5:1	-	1:1	-
MPO CHANNELS 1 to 6	Threshold Kneepoint	-16	OFF	2	OFF	dB
FEATURES	Volume Control	-20	0	2	0	dB
	Power On Delay	0	10	5	0	seconds
	Noise Reduction	OFF	12	6	OFF	dB
	Feedback Canceller	-	-	-	ON	-
	Memory Change Indicator	-	-	-	OFF	-
	Low-Battery Indicator	-	-	-	ON	-


OSPL90 / OSPL50 @ FOG @ IEC 60118-7:2005




OSPL90 / OSPL50 @ FOG @ IEC 60118-0:1993/A1:1994




FREQUENCY RESPONSE @ RTS @ IEC 60118-7:2005




FREQUENCY RESPONSE @ RTS @ IEC 60118-0:1993/A1:1994




INDUCTION COIL SENSITIVITY @ RTS @ IEC 60118-7:2005




FREQUENCY RESPONSE @ VOLUME CONTROL @ IEC 60118-0:1993/A1:1994



INPUT-OUTPUT @ RTS @ IEC 60118-7:2005



INPUT-OUTPUT @ RTS @ IEC 60118-0:1993/A1:1994



--- Damped Ear Hook



Accessories & Spare Parts for Professionals (Fitting)

91118, Brush Cleaner
 53781, Right Programming Cable
 53832, Left Programming Cable
 79345, Battery 13 Flex Adapter
 66183, NOAHLINK Programming Interface (Kernel v. 1.55.03)
 73194, HI-PRO Hearing Instrument Programmer (Firmware 3.00 or higher)
 93769, Neodymium disc magnet
 79430, Membrane Windscreen Bte Artic Grey (10Pcs)
 79431, Membrane Windscreen Bte Beige (10Pcs)
 79432, Membrane Windscreen Bte Artic Grey (10Pcs)
 82876, 6 pack of Microson 13 hearing aid batteries (PR48)
 87304, Membrane Windscreen BTE Black (10Pcs)
 88616, HI-PRO USB Hearing Instrument Programmer (Firmware 4.00 or higher)
 89600, Fitting Software Microson CODA e-STUDIO 6 (6.5.0 or higher)
 94349, P13 Ear Hook Soft rev2 (5pcs)

Accessories & Spare Parts for end user

91118, Brush Cleaner
 94166, User Manual m2 BTE LP2 (ES/EN/IT/PT/FR)
 93769, Neodymium disc magnet
 94583, 6 pack of Microson 13 hearing aid batteries M/Free (PR48)
 88192, Microson Microbox case

HEARING INSTRUMENT CLASSIFICATION IN COMPLIANCE WITH IEC 60601-1 STANDARD

Medical Device Classification

Protection against electric shock	MEDICAL DEVICE WITH INTERNAL ELECTRICAL POWER SOURCE
	Type B Applied Part
	 This symbol indicates that this product adheres to the requirements established for an application component of type B in accordance with IEC 60601-2-66. The surface of the hearing aid is classified as an application component of type B.
Mode of operation	CONTINUOUS OPERATION

Environmental conditions

		Operating	Storage and Transport
	Temperature Limit	0 to 40°C	-20 to 60°C
	Humidity limitation	10 al 95% (Non-condensing)	

Power Supply Electrical Features	m2 BTE
Nominal Operating Voltage	1.4 V
Current Type	Direct current DC
Nominal Current Leakage	0.88 mA
Battery Nomenclature (IEC 60086)	PR48

PRODUCT	REFERENCE	MODEL	GTIN-13
m2 BTE	94115	MICROSON m2 BTE AT BEIGE	8435281312047
	94336	MICROSON m2 BTE AT ARTIC GREY	8435281312030
	94337	MICROSON m2 BTE AT DOLPHIN GREY	8435281312061
	94338	MICROSON m2 BTE AT BLACK	8435281312054

GMDN Code: 34671

Product Data

Measurements performed using a UPL 66 (Rohde & Schwarz) Id 23564 test unit on June 2015 and are subject to changes without prior notice.

