

Earmold for BTE devices


Description

Microson earmold to be connected to any manufacturer and model of BTE hearing aid

Intended Use

The earmold is a custom-made device that is fitted into the natural opening of the patient's ear canal. When used in conjunction with a hearing aid, the device compensates for hearing loss. MICROSON earmolds it from scratch, adapting it to the specific characteristics of the patient's ear canal. It is a reusable device; it does not require fitting. It is designed for adults and children.

The device is not intended for people with intellectual disabilities, due to the high risk of incorrect use. It is a non-mass-produced device, manufactured to the prescription of a qualified healthcare professional, and is intended for use only by a specific patient for the exclusive purpose of satisfying their condition and needs.



Features

Material		Hard	Soft shore 60	Soft shore 40	Soft shore 25	Semisoft
Color		Pink transparent/Transparent				Transparent
Options	Elbow	Available	Not Available			
	Tube	Thin/Standard/Anti-humidity	Thin/Standard/Anti-humidity/with retainer			Thin/Standard/Anti-humidity
	Pull Cord	Available				
	Venting	Ø 1.00 mm - maximum allowed				
	Locking	Canal/Helix/Plus	Canal			



Features

Features		Child		Adult			
Material		Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40	Soft shore 25/Semi soft
Color		Pink transp./Transp.	Not Available	Pink transp./Transp.			Not Available
Options	Elbow	Available		Available	Not available		
	Tube	Thin/Standard/Anti- humidity		Thin/Standard/Anti- humidity	Thin/Standard/Anti-humidity/with retainer		
	Shell	Half / Full		Half / Full			
	Venting	Ø 1mm - Max. allowed		Ø 1 mm – Max. allowed			
	Tinnitus	Available		Available			
	Skeleton	Half / Full		Half/Full			

Product Data

Earmold for BTE devices



Features		Child					Adult				
Material		Hard	Soft shore 60	Soft shore 40	Soft shore 25	Semi soft	Hard	Soft shore 60	Soft shore 40	Soft shore 25	Semi soft
Color		Pink transp./Transp.	Pink transp./Transp. /Blue/Red/ Green/Yellow/ Rose/Fucsia	Pink transp./transp. /Blue/Red/ Yellow/ Turquoise	Not available		Pink transparent/Transparent				Transparent
Options	Elbow	Available	Not available				Available	Not Available			Thin/Standard /Anti-humidity
	Tube	Thin/Standard/A nti-humidity	Thin/Standard/Anti-humidity/with retainer				Thin/Standard /Anti-humidity	Thin/Standard/Anti-humidity/with retainer			
	Shell	Full		Half / Full							
	Venting	Ø 1.00 mm – Max. allowed		Ø 1.00 mm – Max. allowed							

List of BTE earmolds

Part Number	Modelo
EM200077	EARMOULD CANAL SOFT 60
EM200076	EARMOULD CANAL HARD
EM200079	EARMOULD CANAL SOFT 40
EM200181	EARMOULD CANAL SOFT 25
EM200080	EARMOULD CANAL LOCK HARD
EM200081	EARMOULD CANAL LOCK SOFT 60
EM200083	EARMOULD CANAL LOCK SOFT 40
EM200182	EARMOULD CANAL LOCK SOFT 25
EM200157	EARMOULD CANAL SEMI SOFT
EM200161	EARMOULD CANAL LOCK SEMI-SOFT
EM200055	EARMOULD CONCHA HARD
EM200054	EARMOULD SEMI-CONHA HARD
EM200056	EARMOULD CONCHA SOFT 60
EM200084	EARMOULD SEMI-CONHA SOFT 60

Part Number	Modelo
EM200058	EARMOULD CONCHA SOFT 40
EM200045	EARMOULD SEMI-CONHA SOFT 40
EM200183	EARMOULD CONCHA SOFT 25
EM200049	EARMOULD SKELETON SOFT 40
EM200046	EARMOULD SKELETON HARD
EM200050	EARMOULD SEMI-SKELETON HARD
EM200059	TINNITUS EARMOULD HARD
EM200047	EARMOULD SKELETON SOFT 60
EM200158	EARMOULD CONCHA SEMI-SOFT
EM200159	EARMOULD SEMI-CONHA SEMI-SOFT
EM200185	EARMOULD SEMI-CONHA SOFT 25
EM200053	EARMOULD SEMI-SKELETON SOFT 40
EM200051	EARMOULD SEMI-SKELETON SOFT 60
EM200060	TINNITUS EARMOULD SOFT

Product Data

 **WARNING!**

Any modification of the product must be carried out by a Technical Service authorized by the manufacturer.

Discontinued

DS-0100-001-EN
Rev. F 2026-01-26

Earmold for BTE Devices