

Earmold for Open-fit devices




Description


Microson earmold for OPEN FIT hearing aids manufactured by Sivantos, Starkey, GN Resound, Oticon or Phonak.

Intended Use

Earmold is an adapted device that is intended to be placed into the patient ear canal natural orifice. The device combined with hearing aid permits to compensate for impaired hearing. The device is shaped and formed by MICROSON to conform to the patient's ear. This is a reusable device. The device is intended for adults and children. The device is not intended to be used by people with mental disability, due to the high risk of misuse of the device.



Canal Earmold



| Features | | Sivantos | | | Starkey | | Oticon | | Phonak | | GN Resound | |
|----------|-----------|--------------------------|---------------|-----------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|
| Material | | Hard | Soft Shore 60 | Soft shore 40/25/ Semi soft | Hard | Soft shore 60/40/25/ Semi soft |
| Color | | Pink transp./Transp. | | Not available | Pink transp./Transp. | Not available | Pink transp./Transp. | Not available | Pink transp./Transp. | Not available | Pink transp./Transp. | Not available |
| Options | Pull cord | Available | | | Available | | Available | | Available | | Available | |
| | Locking | Available | | | Available | | Available | | Available | | Available | |
| | Venting | Ø 1.00 mm - max. Allowed | | | Ø 1.00 mm - max. Allowed | | Ø 1.00 mm - max. Allowed | | Ø 1.00 mm - max. Allowed | | Ø 1.00 mm - max. Allowed | |
| | k shape | Available | Not avail. | Available | Not available | Not available | Not available | | | | | |

Product Data

Skeleton Earmold



| Features | | Sivantos | | Starkey | | Oticon | | Phonak | | GN Resound | |
|----------|-----------|--------------------------|--------------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|--------------------------|--------------------------------|
| Material | | Hard | Soft shore 60/40/25/ Semi soft |
| Color | | Pink transp./Transp. | | Pink transp./Transp. | Not available |
| Options | Pull cord | Available | |
| | Locking | No disponible | | No disponible | | Available | | Available | | Available | |
| | Venting | Ø 1.00 mm - max. Allowed | |
| | k shape | Available | | Available | Not available | Not available | Not available | | | | |

DS-0100-002-EN
Rev. E 2026-01-26

Earmold for Open-Fit devices

Earmold for Open-fit devices

microson

List of variants

| Part Number | Model |
|-------------|----------------------------|
| 108894 | Molde Open Duro Sivantos |
| 108895 | Molde Open Duro Starkey |
| 112240 | Molde Open Blando Sivantos |
| 108892 | Molde Open Duro GN |
| 108961 | Molde Open Duro Oticon |
| 108893 | Molde Open Duro Phonak |

Product Data

 **WARNING!**

Any modification of the product must be carried out by a Technical Service authorized by the manufacturer.

Dicada

DS-0100-002-EN
Rev. E 2026-01-26

Earmold for Open-Fit devices