

Earmold for Open-fit devices

microson


Description

Microson earmold for OPEN FIT hearing aids.

Intended Use

The earmold is a custom-made device that is fitted into the natural opening of the patient's ear canal. When used in conjunction with a hearing aid, the device compensates for hearing loss. MICROSON earmolds it from scratch, adapting it to the specific characteristics of the patient's ear canal. It is a reusable device; it does not require fitting. It is designed for adults and children.

The device is not intended for people with intellectual disabilities, due to the high risk of incorrect use. It is a non-mass-produced device, manufactured to the prescription of a qualified healthcare professional, and is intended for use only by a specific patient for the exclusive purpose of satisfying their condition and needs..



Canal Earmold



Features		Sivantos			Starkey		Oticon		Phonak		GN Resound	
Material		Hard	Soft Shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft
Color		Pink transp./Transp.		Not available	Pink transp./Transp.	Not available	Pink transp./Transp.	Not available	Pink transp./Transp.	Not available	Pink transp./Transp.	Not available
Options	Pull cord	Available			Available		Available		Available		Available	
	Locking	Available			Available		Available		Available		Available	
	Venting	Ø 1.00 mm - max. Allowed			Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed	
	k shape	Available	Not avail.		Available		Not available		Not available		Not available	

Skeleton Earmold



Features		Sivantos		Starkey		Oticon		Phonak		GN Resound	
Material		Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft	Hard	Soft shore 60/40/25/ Semi soft
Color		Pink transp./Transp.		Pink transp./Transp.	Not available	Pink transp./Transp.	Not available	Pink transp./Transp.	Not available	Pink transp./Transp.	Not available
Options	Pull cord	Available		Available		Available		Available		Available	
	Locking	No disponible		No disponible		Available		Available		Available	
	Venting	Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed		Ø 1.00 mm - max. Allowed	
	k shape	Available		Available		Not available		Not available		Not available	

Product Data

Earmold for Open-fit devices

List of Open-Fit earmolds

Part Number	Modelo
EM200074	EARMOULD OPEN HARD SIVANTOS
EM200075	EARMOULD OPEN HARD STARKEY
EM200163	EARMOULD OPEN SOFT SIVANTOS 3.0
EM200072	EARMOULD OPEN HARD GN
EM200085	EARMOULD OPEN HARD OTICON
EM200073	EARMOULD OPEN HARD PHONAK
EM100686	EMOLDE OPEN DURO WSA

Product Data

 WARNING!

Any modification of the product must be carried out by a Technical Service authorized by the manufacturer.



DS-0100-002-EN
Rev. F 2026-08-13

Earmold for Open-Fit devices