

Earmold for RIC devices




Description

Microson earmold for RIC hearing aids manufactured by Sivantos, Starkey, GN Resound, Oticon or Phonak.

Intended Use

Earmold is an adapted device that is intended to be placed into the patient ear canal natural orifice. The device combined with hearing aid permits to compensate for impaired hearing. The device is shaped and formed by MICROSON to conform to the patient's ear. This is a reusable device. The device is intended for adults and children. The device is not intended to be used by people with mental disability, due to the high risk of misuse of the device.



Canal Earmold



Features		Starkey			Sivantos			Oticon			Phonak			GN Resound		
Material		Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft
Color		Pink transp./Transp.			Pink transp./Transp.			Pink transp./Transp.			Pink transp./Transp.			Pink transp./Transp.		
Opt.	Receiver	50dB, 60dB			S,M,P			60dB, 85dB			XS,S,M,P,XP, S,M,P			LP,MP,HP		
	Pull cord	Available														
	Locking	Available														
	Venting	Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed		
	Open	Avail.	Not Avail.													

Product Data

Skeleton Earmold



Features		Starkey			Sivantos			Oticon			Phonak			GN Resound		
Material		Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft
Color		Pink transp./Tr ansp.			Pink transp./Tr ansp.			Pink transp./Tr ansp.			Pink transp./Tr ansp.			Pink transp./Tr ansp.		
Opt.	Receiver	50dB, 60dB			S,M,P			60dB, 85dB			XS,S,M,P,XP			LP,MP,HP		
	Pull cord	Available														
	Venting	Ø 1 mm - Max. Allowed			Ø 1 mm - Max. Allowed			Ø 1 mm - Max. Allowed			Ø 1 mm - Max. Allowed			Ø 1 mm - Max. Allowed		
	Open	Available														

DS-0100-003-EN
Rev. E 2026-01-26

Earmold for RIC devices

List of variants

Part Number	Model
108885	Molde RIC Duro Sivantos
108882	Molde RIC Duro GN
108889	Molde RIC Duro Oticon
110764	Molde RIC Duro Phonak
112241	Molde RIC Pestaña Duro Sivantos
112245	Molde RIC Pestaña Duro GN
112248	Molde RIC Pestaña Duro Oticon
112249	Molde RIC Pestaña Duro Phonak
108887	Molde RIC Duro Starkey
108888	Molde RIC Blando Starkey
108886	Molde RIC Blando Sivantos
108883	Molde RIC Blando GN
108890	Molde RIC Blando Oticon
112332	Molde RIC Blando Phonak
112243	Molde RIC Pestaña Duro Starkey
112244	Molde RIC Pestaña Blando Starkey
112242	Molde RIC Pestaña Blando Sivantos
112246	Molde RIC Pestaña Blando GN
112247	Molde RIC Pestaña Blando Oticon
112333	Molde RIC Pestaña Blando Phonak

 **WARNING!**

Any modification of the product must be carried out by a Technical Service authorized by the manufacturer.

Diaz