

Earmold for RIC devices

microson


Description

Microson earmold for RIC hearing aids.

Intended Use

The earmold is a custom-made device that is fitted into the natural opening of the patient's ear canal. When used in conjunction with a hearing aid, the device compensates for hearing loss. MICROSON earmolds it from scratch, adapting it to the specific characteristics of the patient's ear canal. It is a reusable device; it does not require fitting. It is designed for adults and children.

The device is not intended for people with intellectual disabilities, due to the high risk of incorrect use. It is a non-mass-produced device, manufactured to the prescription of a qualified healthcare professional, and is intended for use only by a specific patient for the exclusive purpose of satisfying their condition and needs.



Canal Earmold



example image

Features		Starkey			Sivantos			Oticon			Phonak			GN Resound		
Material		Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft
Color		Pink transp./Transp.		Not available	Pink transp./Transp.		Not available	Pink transp./Transp.		Not available	Pink transp./Transp.		Not available	Pink transp./Transp.		Not available
Opt.	Receiver	50dB, 60dB			S,M,P			60dB, 85dB			XS,S,M,P, XP S,M,P					
	Pull cord	Available			Available			Available			Available					
	Locking	Available			Available			Available			Available					
	Venting	Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed			Ø 1.00 mm – Max. Allowed					
	Open	Avail.	Not Avail.		Avail.			Not Avail.			Avail.			Not Avail.		



example image

Skeleton Earmold

Features		Starkey			Sivantos			Oticon			Phonak			GN Resound		
Material		Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft	Hard	Soft shore 60	Soft shore 40/25/ Semi soft
Color		Pink transp./Tr ansp.		Not available	Pink transp./Tr ansp.		Not available	Pink transp./Tr ansp.		Not available	Pink transp./Tr ansp.		Not available	Pink transp./Tr ansp.		Not available
Opt.	Receiver	50dB, 60dB			S,M,P			60dB, 85dB			XS,S,M,P,XP			LP,MP,HP		
	Pull cord	Available			Available			Available			Available			Available		
	Venting	ø 1 mm - Max. Allowed			ø 1 mm - Max. Allowed			ø 1 mm - Max. Allowed			ø 1 mm - Max. Allowed			ø 1 mm - Max. Allowed		
	Open	Available			Available			Available			Available			Available		

Product Data

List of RIC earmolds

Part Number	Modelo	Part Number	Modelo
EM200065	EARMOULD RIC HARD SIVANTOS	EM200063	EARMOULD RIC SOFT GN
EM200062	EARMOULD RIC HARD GN	EM200070	EARMOULD RIC SOFT OTICON
EM200069	EARMOULD RIC HARD OTICON	EM200176	EARMOULD RIC SOFT PHONAK
EM100579	EARMOULD RIC HARD PHONAK	EM200168	EARMOULD RIC CANAL LOCK HARD STARKEY
EM200164	EARMOULD RIC CANAL LOCK HARD SIVANTOS	EM200166	RIC STARKEY SOFT TAB
EM200170	EARMOULD RIC CANAL LOCK HARD GN	EM200165	RIC SIVANTOS SOFT TAB
EM200172	EARMOULD RIC CANAL LOCK HARD OTICON	EM200171	EARMOULD RIC CANAL LOCK SOFT GN
EM200174	EARMOULD RIC CANAL LOCK HARD PHONAK	EM200173	EARMOULD RIC CANAL LOCK SOFT OTICON
EM100595	EARMOULD RIC CANAL LOCK HARD WIDEX	EM200177	EARMOULD RIC CANAL LOCK SOFT PHONAK
EM100529	EARMOULD RIC SKELETON HARD GN	EM100600	EARMOULD RIC CANAL LOCK SOFT WIDEX
EM100523	EARMOULD RIC SKELETON HARD OTICON	EM100526	EARMOULD RIC SKELETON HARD PHONAK
EM100521	EARMOULD RIC SKELETON HARD SIVANTOS	EM100520	EARMOULD RIC SKELETON SOFT GN
EM100527	EARMOULD RIC SKELETON HARD STARKEY	EM100524	EARMOULD RIC SKELETON SOFT OTICON
EM100601	EARMOULD RIC SKELETON HARD WIDEX	EM100522	EARMOULD RIC SKELETON SOFT PHONAK
EM100599	EARMOULD RIC HARD WIDEX	EM100528	EARMOULD RIC SKELETON SOFT SIVANTOS
EM200067	EARMOULD RIC HARD STARKEY	EM100525	EARMOULD RIC SKELETON SOFT STARKEY
EM200068	EARMOULD RIC SOFT STARKEY	EM100594	EARMOULD RIC SKELETON SOFT WIDEX
EM200066	EARMOULD RIC SOFT SIVANTOS	EM100596	EARMOULD RIC SOFT WIDEX

Product Data

 WARNING!

Any modification of the product must be carried out by a Technical Service authorized by the manufacturer.

Dispositiu

DS-0100-003-EN
Rev. F 2026-08-13

Earmold for RIC devices